







ACE® Bearing Combinations

ACE® H-A.C. CUP		BIOLOX® <i>delta</i> CERAMIC	XLPE LINER	XLPE 10° HOODED LINER	XLPE 20° HOODED LINER		DUAL MOBILITY SLEEVE	DUAL MOBILITY XLPE BEARING
36 (BORE)  001-44-00 001-46-00	44/46	 000-44-28	 002-44-28	 003-44-28				
40  001-46-01 001-48-00 001-50-00	46+/48/50	 000-46-32 000-46-28	 002-46-32 002-46-28	 003-46-32 003-46-28	 005-46-32 005-46-28		006-46-36	 007-36-22
44  001-50-01 001-52-00 001-54-00	50+/52/54	 000-50-36 000-50-32	 002-50-36 002-50-32 002-50-28	 003-50-36 003-50-32 003-50-28	 005-50-36 005-50-32 005-50-28		006-50-40	 007-40-22
48  001-54-01 001-56-00 001-58-00	54+/56/58	 000-54-40 000-54-36 000-54-32	 002-54-36 002-54-32 002-54-28	 003-54-36 003-54-32 003-54-28	 005-54-36 005-54-32 005-54-28		006-54-44	 007-44-22 007-44-28
52  001-60-00 001-62-00 001-64-00	60/62/64	 000-60-40 000-60-36 000-60-32	 002-60-36 002-60-32 002-60-28	 003-60-36 003-60-32 003-60-28	 005-60-36 005-60-32 005-60-28		006-60-48	 007-48-22 007-48-28
TRIAL LINERS BLUE = BIOLOX® <i>delta</i> CERAMIC LINER & HIGHLY CROSS-LINKED UHMWPE LINER GREEN = 10° HOODED XLPE LINER ORANGE = 20° HOODED XLPE LINER GREY = DUAL MOBILITY SLEEVE								
XX-XX-XX = Product Code xx-xx-01 = Optimised for large bearings/ Blanking Plugs cannot be used XX-XX-XX = Outside diameter XX-XX-XX = Bearing size + = Large Bore  ACE® LOW PROFILE BONESCREW 009-00- 12, 15, 20, 25, 30, 35, 40, 45, 50, 55, 60								

NB: The use of a skirted femoral head is not recommended for use with ACE® dual mobility